

## GL890 Blue Nitrile™

GL890

Blue nitrile powder free examination glove



### Sectors:

- Automotive • Engineering and Manufacturing
- Food Manufacturing and Preparation
- Janitorial • Medical, Vet and Dental

### Features

**Sensitive Skin:** No latex proteins - eliminates protein sensitisation

**Strength:** Synthetic rubber offers better puncture and abrasion resistance than natural rubber or vinyl

**Comfort:** Soft elastic material feels cooler and is more comfortable and flexible to wear

**Powder Free:** Reduces risk of product contamination

**Grip:** Textured fingertips for improved grip

**Food Contact:** Tested in accordance with European standard EN1186 making it suitable for food contact

Code	GL890
Colour	Blue
Sizes	XS-XXL



Inners dimensions (mm)  
235 x 117 x 78 mm

Outers dimensions (mm)  
418 x 258 x 251mm

Size	Inner Barcode	Outer Barcode
XXS/6	5010699528322	05010699839336
S/7	5010699528339	05010699839343
M/8	5010699528346	05010699839350
L/9	5010699528353	05010699839367
XL/10	5010699536570	05010699839374
XXL/11	5010699536570	05010699849809

### Standards

Contains accelerators which may cause allergic reactions

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with the Medical Devices Regulation (EU) 2017/745 Class 1

Tested to European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Product Food Approved to EC1935/2004

EU Type Examination carried out by SATRA

Tested to European Standards EN420:2003+A1:2009 & EN ISO 374-1:2016 & EN ISO 374-5:2016

### Packing

100 gloves per dispenser  
10 dispensers per case  
50 cases per pallet



MDR CLASS I  
PPE CAT III  
2777



K



VIRUS



For further technical information, please contact Customer Services on 033 3320 8550 or the Technical team on 033 3320 8555

